


CATALOG NO.	* 3005-42/-70	
SIZE (NOM VOL. (mL)	500	
VOL. AT NECK CAP. (mL)	670	
BOTTLE	PETG	
CLOSURE	WHITE PP	
CLOSURE LINER	LDPE FOAM	
NECK FINISH	38-430	
THREAD (T.P.I.)	6	
NOM. WEIGHT W/ CLOSURE (gms)	85	
	INCH	CM
A (NECK I.D.)	1.11	2.8
B (HT. W/O CLOSURE)	6.80	17.2
C (HT. W/ CLOSURE)	6.94	17.6
D (SQUARE)	3.0	7.7
W.T. (NOM. WALL THK)020	.05

* 3005-42 IS LAB-PACK, STERILE
3005-70 IS BULK-PACK, STERILE

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Pancosmo Creek Drive
Rochester, NY 14625-2385, USA

NAME MC² PVC
CR# PVC

APP'D 1/8/06

REV 0106

PSS 3005-XX

BIOTAINER BOTTLE, PETG, 500mL

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.