



Bel-Art · H-B Instrument · Wilmad-LabGlass

661 Route 23 South, Wayne, New Jersey 07470

Tel: 973-694-0500 * Fax: 973-694-7199

Internet Address: <http://www.bel-art.com>

CERTIFICATE OF ANALYSIS

This is to verify that all items noted below are in conformance with the Purchase Order Contract, drawings, specifications, and other applicable documentation; that all process certifications, material certifications, are on file at this or another Bel-Art facility and are available for review by the customer and or Government Representatives upon request.

Customer: Spectrum Lab Products

P.O. #: 1233854

Order #: B441362

Product #: 37260-0000 **Mfg. Date:** 9/8/2017 **Lot #:** J223897

Description: Liquid Nitrogen Cooled Mortar

Material: Polyethylene, EPS Foam, 304 Stainless steel

Dimensions: OD 6 1/2" (165mm) Diameter, 4 1/2" (114mm) High
ID of U shape bowl, 4" (101mm), 1.9" (48mm) Radius
Removable Handle 4" long

Other: Pestle Not Included
Any Pestle that is appropriate for use with a 275ml mortar,
Example: Coors Porcelain Pestle 60320

Manufactured: Bel-Art Products: Wayne N.J.; U.S.A.

Reviewed and Approved By:

Signature: *Michael E Grzelak*

Printed Name: Michael E Grzelak

Title: Document Control Coordinator/Quality Analyst

Date: 9/13/2017

SUBSIDIARIES: MADDAK, INC * APPLIED COATING, INC
ADDITIONAL PLANTS IN POCOMOKE, MD * KANSAS CITY, MISSOURI