

TECHNICAL DATA SHEET

Part Nbr	088783
Designation	Cooling case Cryo-safe cooler



Description :

- ✓ PC refrigerant case filled with a non-toxic gel
- ✓ For 12 tubes of 15 ml
- ✓ Holds tubes at a temperature below 0°C for approximately 3 hours
- ✓ Rubber feet prevent sliding on the lab bench
- ✓ Stackable
- ✓ Color : orange – transparent lid
- ✓ Weight : 1,6 kg
- ✓ Dimensions : 19,7 x 14 x 19 cm

Sales unit :

1 Cooling case Cryo-safe cooler
