



**AQUET DETERGENT
CONCENTRATE FOR GLASSWARE AND
PLASTICS; 20 ML POUCHES (PACK OF
20)**

Phosphate-Free Liquid Detergent with Neutral pH

Catalog No. : F17094-0020
UPC No. : 191634006404

Safe for glassware, plastics and surfaces, this phosphate-free, nonionic surfactant will not etch or mar precision glassware. It rinses quickly and dries without wiping.

- Convenient, pre-portioned pouches make one gallon of 1% solution
- 3.8 liter (1 gallon) bottle comes with 473ml (16oz) working solution dispenser bottle.
- 20 liter (5 gallon) Cubitainer™ for frequent users
- Not for use in automatic dishwashers

Catalog No.:	Size	Pack	Quantity / Pack
F17094-0020	20 ml pouch makes one gallon	20 pouches	20

DOMINIQUE DUTSCHER SAS