

Cryo Tube



Features of Cryo Tube:

- Storage of samples in freezers or similar low-temperature equipment
- External thread without silicone gasket
- Secure and tight closure of the tube
- Self-standing
- Single-handed operation with the use of the cryogenic rack No. 99016 with its star-shaped locking system in the grooves
- Color-coding inserts
- Convenient, puncture and tear resistant re-closable bag (picture 4, following page)

During storage build-up of stress on the product can happen. The TPP design as well as the high grade material ensures safe handling during and after the thawing process.

- 1 Smooth internal surface eliminates undue loss of specimen material during removal.
- 2 Lid can be removed single handed with ease, even after multiple usages.






- 1 The yellow cryo rack #99016.
- 2 A single-handed operation is possible using the tubes with the cryogenic rack # 99016 with its star-shaped locking system in the grooves.

Do not use cryo tubes for storage in the liquid phase but in the vapor-phase above the level of the liquid-phase nitrogen. During storage in liquid nitrogen TPP recommends use of additional wrapping such as welded tubes. Always use appropriate safety equipment and procedures.

Type	Product-No.	Volume Skala max. ml	Volume recom. ml	Dimensions Ø x h mm	Colours	Material	Qty / Bag	Qty / Case
	89012	0.9	0.7	12 x 37	-	PP	100	800
	89020	1.5	1.3	12 x 48	-	PP	100	800
	89040	3.5	3.3	12 x 75	-	PP	100	400
	89050	4.5	4.3	12 x 90	-	PP	100	400
	99020	-	-	11 x 1	6 colours	PP	6 x 100	7800
	89801	-	-	11 x 1	white	PP	600	3600
	89802	-	-	11 x 1	blue	PP	600	3600
	89803	-	-	11 x 1	yellow	PP	600	3600
	89804	-	-	11 x 1	green	PP	600	3600
	89805	-	-	11 x 1	pink	PP	600	3600
	89806	-	-	11 x 1	red	PP	600	3600

Type	Product-No.	Fitting Vials	Capacity Qty	Dimensions l x w x h mm	Material	Qty / Bag	Qty / Case
	99014	89040 – 89050	81	133 x 133 x 95	PP	1	10
	99015	89012 – 89020	81	133 x 133 x 46	PP	1	20
	99016	89012 – 89050	40	100 x 200 x 25	PC	1	12



- 3 Cryo tubes with volume scale and writing area.
- 4 Convenient re-closeable bag with zipper.