

Product Specification

Product Code	BDA212P
Product Description	AZLON BEAKER GRADUATED 25ML PP
Material particular	Polypropylene
Capacity ml	25
Colour	Natural
Pack quantity	10

Dimensions and tolerances

Height mm	40
Diameter mm	42
Graduation ml	2

Standards

Conforms to	ISO 7056
-------------	----------

Product Compliance

Reach	Yes
-------	-----

SciLabware Limited hereby confirms that the product specified above is inspected according to all relevant specifications and procedures laid out in our Quality System, in accordance with our registration to BS EN ISO 9001:2008

Edel Tuite

Azlon Quality Manager

Date issued 4/12/2016