


PRODUCT SPECIFICATION

Product Code:	3151/14
Product Description:	100ml Volumetric Class A Flask
Brand:	PYREX
Type of Glass:	Borosilicate glass Type 3.3, which conforms to ISO3585 European and US Pharmacopoeia Type 1 glass
Product Specification:	Finely calibrated to USP standard Individual serial number, batch code and date of manufacture laser etched Downloadable batch certificate can be obtained from our website www.scilabware.com Supplied with polyethylene stopper Conformity 'H' certified
	Tolerance ±0.08ml Approx. height 170mm Stopper size 14/23
Product Standard(s):	The product is manufactured to comply with ISO1042 and DIN 12664
Pack Quantity:	5 pieces per pack

SciLabware Limited hereby confirms that the product specified above is inspected according to all relevant specifications and procedures laid out in our Quality System, in accordance with our registration to BS EN ISO9001 and BS EN ISO13485.

Signed on behalf of SciLabware Limited by:



Signed: _____

K. Skupniewicz
Quality Engineer

Date: 13 Feb 2018