

Product Specification

Catalogue No:	3151/16
Product Description:	200ml Volumetric Flask. Class 'A' USP/ISO/DIN
Brand:	PYREX
Type of Glass:	Borosilicate Glass; Type 3.3, which complies with ISO3585 European and United States Pharmacopoeia Glass Type 1
Product Markings:	Product branding and marking spot, volumetric capacity, tolerance and related information.
Batch Reference code:	For traceability the flask is laser etched with an individual serial number and date on the base. A certificate of performance for the batch can be downloaded from our website - www.scilabware.com
Neck/ Stopper:	Supplied with a 14/23 White Polyethylene Stopper
Specification:	Volumetric tolerance = +/- 0.10ml Overall height = 210mm approx. External diameter, flask body = 76mm approx. Socket size = 14/23
Standard:	Conforms to ISO1042 and DIN12664
Pack Quantity:	Quantity 5 pieces per pack

SciLabware Limited confirms that the above mentioned product has been manufactured to meet all relevant standards as detailed by our Quality Management System.



R J Turner
(Quality Manager)