

PRODUCT SPECIFICATION

Product Code:	CRC209
Product Description:	100ml Measuring Cylinder - Class 'B' with blue graduations
Brand:	MBL
Glass Type:	Borosilicate Glass Type 3.3 which complies with ISO3585 European and United States Pharmacopoeia Glass Type 1
Specification:	Graduation 1.0ml Tolerance ± 1.0 ml Max height 260mm
Product Standard(s):	The product is manufactured to comply with ISO4788
Pack Quantity:	2 pieces per pack

SciLabware Limited hereby confirms that the product specified above is inspected according to all relevant specifications and procedures laid out in our Quality System, in accordance with our registration to BS EN ISO9001:2008.

Signed on behalf of SciLabware Limited by:



Signed: _____

K. Skupniewicz
Quality Engineer

Date: 23 Jan 2018