

## Product Specification

Product Code	CPB0500P
Product Description	AZLON MEASURING CYLINDER, CLASS B, BLUE PRINTED GRADUATIONS, PP, 500 ML
Material particular	Polypropylene
Capacity ml	500
Colour	Ultra Clear PP
Print Colour	Blue
Autoclavable	Yes, see <a href="http://www.scilabware.com">www.scilabware.com</a>
Pack quantity	1

## Dimensions and tolerances

Height mm	360
Diameter mm	55
Graduation ml	50.0

## Standards

Conforms to	ISO 6706
-------------	----------

## Product Compliance

TSE/BSE Compliant	Yes
Reach	Yes
Latex	Yes

SciLabware Limited hereby confirms that the product specified above is inspected according to all relevant specifications and procedures laid out in our Quality System, in accordance with our registration to BS EN ISO 9001:2008

Edel Tuite

Azlon Quality Manager

Date issued 1/19/2017