

## Description

- Class A volumetric flasks with blue markings 250 ml

## Part number

- 090236

Pieces by sales : 2



---

## Description

- Printed with blue enamel high contrast markings, these Class A volumetric flasks comply with the requirements of ISO and DIN standards.
- Complies with ISO 1042 and DIN 12664
- Finely calibrated to Class A tolerances
- Each flask is laser etched with an individual serial number batch code and date of manufacture
- Manufactured from chemically resistant borosilicate glass
- Conformity 'H' certified
- Inscriptions in blue enamel
- Supplied with polyethylene stopper
- Batch certificate available to download online
- Approx. height 220 mm
- Capacity 250 ml
- Stopper size 14/23