

► Thermo Scientific Sterilin BenchGuard



Sterilin BenchGuard is highly absorbent paper ideal for protecting benches and other surfaces against liquid spills.

details

- Choice of absorption capacity — 400mL/m² or 800mL/m²
- One side plastic coated to prevent soak through
- Available in roll or sheet form
- Rolls supplied in easy-to-use dispenser packs

Ordering Information: Not available in all geographic regions

Sterilin BenchGuard

Cat. No.	Description	Absorption	Length x Width	Pack Type
BG50	BenchGuard	400mL/m ²	50.0 × 0.49m	1 x Roll
BG60/100	BenchGuard	400mL/m ²	0.60 × 0.49m	100 x Sheets
BG50E	BenchGuard Extra	800mL/m ²	50.0 × 0.49m	1 x Roll
BG60E	BenchGuard Extra	800mL/m ²	0.60 × 0.49m	50 x Sheets

DOMINIQUE DUTSCHER SAS