

TECHNICAL DATA SHEET

Part Number	090478
Description	Pyrex® Culture tubes



Benefits :

- ✓ Lightweight microtube racks are ideal for long-term studies or tube storage.
- ✓ Two-piece racks are constructed of expanded polystyrene (EPS) foam
- ✓ Extra spacing between wells allows easy gripping with fingertips
- ✓ Racks stack securely with nesting features on lid and base
- ✓ Ideal for culture work
- ✓ Manufactured from Pyrex® borosilicate glass
- ✓ High resistance to attack from water reduces leaching of contaminants which can cause pH changes
- ✓ PBT screwcaps with inert PTFE lined rubber disc
- ✓ Cap can withstand repeated dry heat and steam sterilisation cycles.

Technical specifications :

- ✓ Capacity : 11 ml
- ✓ OD x length : 16 x 100 mm
- ✓ Cap size : 15 mm
- ✓ Wall thickness \pm 0.2mm : 1.8

Sales unit :

40 Culture tubes