



SciLabware

QF/7.17 Issue A 10 Nov 2017

SciLabware Limited
Unit 4, Riverside 2, Campbell Road
Stoke-on-Trent, ST4 4RJ, United Kingdom

Tel +44 (0) 1782 444406

Fax +44 (0) 1782 940436

www.scilabware.com
www.wheaton-uk.com

Registered office (as above)
Registered in England and Wales No. 7723574
BS EN ISO 9001: 2008 registered company

PRODUCT SPECIFICATION

| | |
|-------------------------------|---|
| Product Code: | 1636/04MP |
| Product Description: | 14ml Culture Tube |
| Brand: | PYREX |
| Material | |
| Glass Tube: | Clear Borosilicate glass Type 3.3, which conforms to ISO3585 European and US Pharmacopoeia Type 1 glass |
| Screw Cap: | Black Polybutylene terephthalate (PBT) |
| Liner: | PTFE faced ethylene propylene elastomer |
| Product Specification: | |
| Glass Tube: | Tube O.D. 16mm Tube Length 125mm Tube Wall thickness 1.8mm Thread size 15mm |
| Screw Cap: | Cap O.D. 21.4mm Cap Height 19mm |
| Liner: | Liner O.D. 15.2mm Liner Thickness 2mm Max. |
| Pack Quantity: | 40 pieces per pack |
| Manufacturer: | DWK Life Sciences, Germany |

SciLabware Limited hereby confirms that the product specified above is inspected according to all relevant specifications and procedures laid out in our Quality System, in accordance with our registration to BS EN ISO9001 and BS EN ISO13485.

Signed on behalf of SciLabware Limited by:

Signed: M.Ashley

Quality Technician

Date: 05 Jul 2019

Pyrex® is a registered trademark. SciLabware registered user.

Azlon®

MBL®

PYREX®



SVL®

WHEATON®