

## PRODUCT SPECIFICATION

<b>Product Code:</b>	CRC204
<b>Product Description:</b>	25ml Measuring Cylinder - Class 'B' with blue graduations
<b>Brand:</b>	MBL
<b>Glass Type:</b>	Borosilicate Glass Type 3.3 which complies with ISO3585 European and United States Pharmacopoeia Glass Type 1
<b>Specification:</b>	Graduation 0.5ml Tolerance $\pm 0.5$ ml Max height 170mm
<b>Product Standard(s):</b>	The product is manufactured to comply with ISO4788
<b>Pack Quantity:</b>	2 pieces per pack

SciLabware Limited hereby confirms that the product specified above is inspected according to all relevant specifications and procedures laid out in our Quality System, in accordance with our registration to BS EN ISO9001:2008.

Signed on behalf of SciLabware Limited by:



Signed: \_\_\_\_\_

K. Skupniewicz  
Quality Engineer

Date: 23 Jan 2018