

Data sheet

Product: ringcaps® 44,7 µl, sodium heparine

Code-no.: 960 21 44

The above named product conforms with the requirements of the standard ISO 7550.

Glass material: DURAN, borosilicate glass 3.3 (acc. ISO 3585)

Inner diameter: 1,02 ± 0,05 mm

Outer diameter: 1,52 ± 0,05 mm

Tolerances: Accuracy A: ±0,25%
Coefficient of variation: 0,50%

Graduation line: black, at the volume 44,7 µl

Eberstadt, 08.10.2018



.....
Roland Benz/ Quality Manager