

Technical data sheet

Article: haematocrit tubes
Code no.: 910 01 ..
910 02 ..

Material:

Soda-lime glass

Chemical properties:

- Hydrolytic resistance at 98°C: Class 3 acc. ISO 719
- Hydrolytic resistance at 121°C: Class 3 acc. ISO 720
- Acid resistance: Class 1 acc. ISO 1776
- Alkali resistance: Class 2 acc. ISO 695

Physical properties:

- Transformation temperature: 525°C
- Density: 2,50 g/cm³
- Linear expansion coefficient: 9,1 x 10⁻⁶ K⁻¹

Dimensions:

Code no:	length	outer diameter	volume
9100160, 9100260	75 ±0,5 mm	1,3 - 1,4 mm	60 µl
9100170, 9100270	70 ±0,5 mm	1,5 - 1,6 mm	70 µl
9100175, 9100275	75 ±0,5 mm	1,42 - 1,68 mm	75 µl

Hirschmann Laborgeräte
GmbH & Co.KG



Roland Benz
Quality Manager

Eberstadt, 09.02.2015