

Data sheet

Article: Measuring cylinder 10 ml Class A

Code no.: 224 01 60

Material:

DURAN (Borosilicat glass 3.3)

Chemical properties:

- Hydrolytic resistance at 98°C: Class 1 acc. ISO 719
- Hydrolytic resistance at 121°C: Class 1 acc. ISO 720
- Acid resistance: acc. ISO 1776
- Alkali resistance: Class 2 acc. ISO 695

Nominal volume:

10 ml

Subdivision of the scale:

0,2 ml

Tolerance:

± 0,1 ml

Hirschmann Laborgeräte
GmbH & Co.KG



Roland Benz
Quality Manager

Eberstadt, 17.06.2019