

TECHNICAL DATA SHEET

Part Number	090589 to 090598
Description	Bulb pipettes



Benefits :

- ✓ Made from soda-lime glass AR-Glas®
- ✓ 2 marks
- ✓ Waiting time 5 seconds

Certificate :

- ✓ BlauBrand®
- ✓ Class AS
- ✓ DIN EN ISO 648
- ✓ Delivered with batch number and certificate

Technical specifications :

Part number	Nominal volume	Color coded acc ISO 1769	Length	Error limit
090589	0,5 ml	2x noir	300 mm	± 0,005 ml
090590	1 ml	bleu	300 mm	± 0,008 ml
090591	2 ml	orange	330 mm	± 0,01 ml
090592	3 ml	noir	330 mm	± 0,01 ml
090593	5 ml	blanc	400 mm	± 0,015 ml
090594	10 ml	rouge	440 mm	± 0,02 ml
090595	15 ml	vert	510 mm	± 0,03 ml
090596	20 ml	jaune	510 mm	± 0,03 ml
090597	25 ml	bleu	520 mm	± 0,03 ml
090598	50 ml	rouge	540 mm	± 0,05 ml

Sales unit :

1 Bulb pipette