

## Product Specification

<b>Product Code:</b>	1515/08D
<b>Product Description:</b>	Reagent Bottle 500ml - supplied with PBT cap and PTFE faced sealing disc.
<b>Brand:</b>	PYREX
<b>Glass Type:</b>	Borosilicate Glass Type 3.3 which complies with ISO3585 European and United States Pharmacopoeia Glass Type 1
<b>Graduations</b>	White enamel graduations from 100ml up to 400ml in 100ml increments.
<b>Cap:</b>	Polybutylene Terephthalate screw cap with PTFE faced sealing disc
<b>Product Specification:</b>	Capacity – 500ml Bottle height – 175mm Body diameter – 81mm Thread size – SVL30
<b>Standards:</b>	The bottles are manufactured to comply with the requirements of Product Standard ISO4796-1.
<b>Pack Quantity:</b>	10 piece per pack

SciLabware Limited hereby confirms that the product specified above is inspected according to all relevant specifications and procedures laid out in our Quality System, in accordance with our registration to BS EN ISO9001:2008.

Signed on behalf of SciLabware Limited by:



Signed: \_\_\_\_\_  
K. Skupniewicz  
Quality Engineer

Date: 23 Aug 2017