

FICHE TECHNIQUE

Référence	090919
Désignation	Pipette



Description :

- ✓ Graduated cotton plugged pipettes made from high transparency polystyrene
- ✓ Precision: +/- 1 %
- ✓ Endotoxin level less than 0.01 EU/ mL
- ✓ Standard colour code, double graduation
- ✓ Batch number on each box
- ✓ Negative graduations for excess capacity
- ✓ Sterilised by gamma rays, in individual paper/plastic packaging, all plastic packaging or multiple sachets

Unité de vente :

1 Sterile serology pipette
