

CAPP

D E N M A R K

Microbiology

SINGLE CHANNEL PIPETTES

microbiology.capp.dk

CAPP MICROBIOLOGY PIPETTES: IDEAL FOR MICROBIOLOGY TESTS IN THE FOOD INDUSTRY, THIS PIPETTE IS SUPPLIED WITH THREE FIXED VOLUME KNOBS AT 100UL, 900UL AND 1000UL AND HAS A MODIFIED TIP CONE SUITED FOR USE WITH A STERILE EXTENDED LENGTH STRAW, ALLOWING FURTHER REACH AND PIPETTING OF VISCOUS LIQUIDS SUCH AS MILK, CREAM, YOGHURT AND COSMETICS. IT IS AN IDEAL INSTRUMENTS FOR PROTOCOLS INVOLVING DILUTION OF 1:10.

MAIN FEATURES:

- EXTENDED REACH*
- ROTATION FREE VOLUME SELECTION*
- PRE-SET CALIBRATED VOLUME CONTROLLERS*
- FULLY AUTOCLAVABLE*



Ordering Info

CappAero Microbiology pipettes (for use with sterile straws)

CAT. NO.	DESCRIPTION	INACCURACY %	IMPRECISION %
C1000-1BS	CAPPAERO MICROBIOLOGY PIPETTE, WITH 100UL, 900UL AND 1000UL VOLUME KNOBS	3.00 / 1.60	0.80 / 0.30
4170807	PRE-STERILE STRAWS FOR CAPPAERO MICROBIOLOGY PIPETTE, 40X25 PCS.		