

Specification

HEPES for molecular biology

A3724

Physical Description:	Solid
Product Code:	A3724
Product Name:	HEPES for molecular biology
Specifications:	DNases/RNases/Proteases: not detectable Assay: min. 99.5 % pH (1 %; H ₂ O): 4.7 - 6.0 Water: max. 0.25 % Chloride: max. 0.005 % Sulfate: max. 0.02 % As: max. 0.0001 % Fe: max. 0.001 % Pb: max. 0.001 % A (1 cm/0.1 M HPLC water) 260 nm: max. 0.1 280 nm: max. 0.05
WGK:	1
Storage:	RT
Molecular Formula:	C ₈ H ₁₈ N ₂ O ₄ S
M:	238.31 g/mol
CAS:	7365-45-9
EINECS:	230-907-9
CS:	29335995

AppliChem GmbH