

# Water Sterile

## Product specificities

Appearance	Clear colorless liquid
Chloride (Cl)	max. 0.00001%
Resistivity (at manuf.)	18-18.5MohmXcm
Iron (Fe)	max. 0.00001%
Heavy metals (as Pb)	max. 0.000001%
Residue after evaporation	max. 0.0001%w/w
Subs. reducing KMnO <sub>4</sub>	Passes test
Sulfate (SO <sub>4</sub> )	max. 0.00005%
TOC	max. 30ppb
DNase activity	Not detected
RNase activity	Not detected
Protease activity	Not detected

Cat. n° **232187**

Formula : H<sub>2</sub>O

CAS : 7732-18-5

EC : 231-791-2



### Standard pack

23218721	6 x 1l	Glass Bottle 28 (1l)
23218731	1 x 100ml	Glass Bottle (100ml)

### More informations

EU FR DE NL IT