


These Petri dishes are made from DURAN® borosilicate 3.3 glass using a special manufacturing process that permits the uniform distribution of agar and guarantees distortion-free observation.

Typical applications: biological and clinical research, cultivation of microorganism, microscopy.

Cat. No.	d ₁ (OD) (mm)	d ₂ (OD) (mm)	h (mm)	h ₁ (mm)	Pack Unit
217554101	60	54	22	20	10
217554307	80	74	22	20	10
217554607	100	94	17	15	10
217554804	100	94	22	20	10
217555106	120	114	22	20	10
217555303	150	143	32	30	10

DUROPLAN® Petri Dish

lid and base have planar surfaces, and free from bubbles and streaks



Cat. No.	d (OD) (mm)	h (mm)	Pack Unit
237553903	40	12	10
237554008	60	15	10
237554205	80	15	10
1184071	90	15	10
237554505	100	10	10
237554608	100	15	10
237554805	100	20	10
237555107	120	20	10
237555201	150	25	10
237555604	180	30	10
237555904	200	30	10
237556103	200	45	10

STERIPLAN® Petri Dish

from Soda-lime Glass

