

TECHNICAL DATA SHEET

Part Number	096035
Description	DURAN® Bromobutyl Rubber Stopper, for GL 45



Benefits :

- ✓ Bromobutyl rubber elastomer conforms with EP 3.2.9 Type I closure for chemical and penetration tests
- ✓ Gas tight seal provides controlled environment in the bottle
- ✓ A Type III Packaging Material DMF (35811) for these products has been filed with the FDA
- ✓ Elastic material for effective sealing and reliable self-sealing after needle piercing (21-gauge "21G" or smaller)
- ✓ Multiple punctures providing easy access to the contents with a syringe
- ✓ Optimised to fit DURAN® GL 25 or GL 45 bottles
- ✓ Suitable for autoclaving / steam sterilising at 121 °C for 20 minutes
- ✓ Wide usable temperature range: -50 to +121 °C

Technical specifications :

- ✓ Color Grey
Diameter d : 41 mm. H : 21 mm. B : 41 mm
- ✓ Thread GL 45
- ✓ Material BIIR
- ✓ Width : 41 mm. Height : 21 mm
- ✓ Stopper Style : straight

Sales unit :

10 Rubber stopper