



CLICK
AND
CARRY

THE NEW
DURAN[®]
GL 45 BOTTLE
CARRYING SYSTEM

DURAN® GL 45 BOTTLE CARRYING SYSTEM

The new modular Bottle Carrying System fits standard laboratory bottles that have a GL 45 neck finish and is compatible with four bottle sizes 2, 5, 10 and 20 litres. It is ergonomically designed, and offers health & safety benefits to the user and allows easy and safe transport of bottles around the laboratory or production facility. It makes pouring liquids from the larger bottles more comfortable and safer. All materials are fully autoclavable.

Carrying Handle

- Carrying handle is suitable for bare or gloved hands
- Handle strap is flexible enough to allow pouring

Tightening Strap

- With indicator arrow
- End modified to prevent unthreading

Pouring Handle

- Ergonomically shaped
- Anti-slip grip



Silicone Support Base

- Structural support for base of larger 5, 10 and 20 litre bottles
- Protection of glass bottle base from impact during lifting or handling
- Spacer to protect glass from contact with Metal Dolly
- Moulded-in holes minimize liquid traps and allow maximum steam penetration during autoclaving



GL 45 Bottle Collar


- The collar is securely attached by pushing down over the GL 45 thread with 4 teeth that lock under transfer ring.
- The collar is only removable by cutting – so attachment is 100% secure in use. A cut line indicator for removal of collar when required.

Horizontal Straps

- Durable stitched construction
- Made from "anti-slip" material



Lot Specific Retrace Code
for traceability and certification



Optional Extra: Metal Dolly

- Stainless steel dolly for 10 or 20 litre bottles (2 sizes only)
- Fully compatible with carrying system
- Easy to clean, crevice-free construction

FOR EASY AND SAFE CARRYING & POURING



Easy and safe carrying using handle

- Handle attaches at the neck and fully supports the base
- For 2, 5, 10 and 20 litre bottles



Safe and controlled pouring using handle

- Safer pouring of liquids using a second lower "Pouring Handle" located on the side strap
- Available on the 5, 10 and 20 litre carriers only

Ordering information

Description	Pack Unit	Cat. No.
DURAN® GL 45 Bottle Carrying System for the 2 Litre DURAN® Bottle (Blue)	1	292836307
DURAN® GL 45 Bottle Carrying System for the 2 Litre DURAN® Bottle (Purple)	1	292826306
DURAN® GL 45 Bottle Carrying System for the 2 Litre DURAN® Bottle (Red)	1	292816305
DURAN® GL 45 Bottle Carrying System for the 2 Litre DURAN® Bottle (Yellow)	1	292806304
DURAN® GL 45 Bottle Carrying System for the 5 Litre DURAN® Bottle (Yellow)	1	292807309
DURAN® GL 45 Bottle Carrying System for the 10 Litre DURAN® Bottle (Yellow)	1	292808605
DURAN® GL 45 Bottle Carrying System for the 20 Litre DURAN® Bottle (Yellow)	1	292809104
DURAN® Metal Dolly for 10 L DURAN® Bottle (stainless steel) available July 2018	1	299018607
DURAN® Metal Dolly for 20 L DURAN® Bottle (stainless steel) available July 2018	1	299019106

* Bottles supplied separately

Fast and secure tightening

of the polypropylene straps with two side strap adjusters for the larger 5, 10 and 20 litre bottles

4 sizes available:
2 L / 5 L / 10 L / 20 L



2 L bottle carrier in 4 colors for color coding





**DURAN
WHEATON
KIMBLE**

Excellence in your hands

DURAN
WHEATON
KIMBLE
DWK Life Sciences GmbH
Dill-Seibert-Str. 21
97077 Wertheim/Main
Germany
www.DWK-LifeSciences.com

RETRACE CODE
DURAN
WHEATON
KIMBLE

DWK Life Sciences GmbH
Hattenbergstraße 10
55122 Mainz, Germany

Phone: +49 6131-1445 4131
Fax: +49 6131-1445 4016
sales@DWK-LifeSciences.com
www.DWK-LifeSciences.com

6038196 EN 1.0 01/2018 Printed in Germany