

TECHNICAL DATA SHEET

Part number	097819
Description	ClearLine reagent reservoirs



Description :

- ✓ White medical grade polystyrene
- ✓ Certified free from RNase, DNase and human DNA, non-pyrogenic
- ✓ Each corner has a pouring spout
- ✓ Graduations on both inner walls (precision +/-5 %)
- ✓ V bottom for maximum sample recovery
- ✓ Non-sterile
- ✓ Manufactured in a clean environment, packaged in a cleanroom
- ✓ Dimensions interior : 117 x 49 x 23 mm
- ✓ Dimensions exterior : 145 x 82 x 27 mm

Pieces by sales :

100 Reservoirs
