

## LCD DIGITAL INDOOR/OUTDOOR THERMOMETER AND INDOOR HYGROMETER



### Réf 9222AT

- Temperature range (In / Out) : -50° to +70°C
- Humidity range (In) : 20% to 100%RH
- Resolution 0,1°, 1%
- Accuracy  $\pm 1^{\circ}\text{C}$ ,  $\pm 3\%$
- Max/Min record
- Stainless probe / Wire 3 m
- Stand / Wall mount
- D : 80x110x20mm / 200g
- Batterie AAA

### Réf 9222AT/111LM – Digital thermo/hygromètre 9222AT delivered with thermal retarder

- Plastic bottle / Alcoholic solution 5 cl
- Rubber seal / Clamping plug
- Ideal to avoid large variations in measured temperature



+

