

Room temperature data logger only

Reference 0176BL

Additional temperature recorder -40 ° C to + 70 ° C, allowing wireless transmission of data to a hub receiver, compatible with our following Recorder + Hub packs:

- Temperature recorder + Hub 176BL / PACK2 pack
- Temperature recorder pack and hygrometry + Hub 177BL / PACK2
- Pack temperature recorder with wired probe "cold chain" + Hub 178BL / PACK2
- Pack temperature recorder wired probe "hot chain" + Hub 179BL / PACK2



- 185,000 memories / Encrypted data
- Temp. : -40.0 ° C to + 70.0 ° C / Accuracy: 0 ° C $\leq \pm 0.2 \text{ ° C}$ <math>< 30 \text{ ° C}</math> ~ $\pm 0.5 \text{ ° C}</math>$
- Current values, historical dated
- Unique serial number: Identification number / QR code
- Sampling intervals: 10 minutes / Up to 3 years of registration
- Waterproof IP67
- Complies with HACCP / EN12830
- D: 85x55x5 mm / 30 g / Wall support
- Battery: Lithium / Autonomy 3 years

