



PRODUCT DATA SHEET

COWIE TECHNOLOGY GROUP LTD
RIDGEWAY, COULBY NEWHAM
MIDDLESBROUGH, CLEVELAND
TS8 0TQ, ENGLAND

TEL: +44 (0)1642 599190
FAX: +44 (0)1642 596810

enquiries@cowie.com
www.cowie.com

Item:	PTFE Bottle
Product Code/Ref:	015.1000.2
No. Pages:	1

Material: PTFE – conformances below.

Capacity: 1litre

Nominal Dimensions: Height 185mm x OD. 100mm x Mouth ID. 57mm.

Max. use Temp. PTFE Withstands 280°C. but the vessel MUST NOT be heated with volatile contents with the lid in place.

Sterilisation: Chemical or thermal but NOT gamma radiation.

Chemical resistance: Virtually total.

GRC: 14.03.17

PTFE used in the manufacture of this product is VIRGIN PTFE GRADE TF1620 manufactured by 3M® DYNEON®

3M® DYNEON® TF1620 – COMPLIANCE STATUS

FDA – Directive FDA 21 CFR 177.1550
USP Class VI – USP 27, NF 22, 2004 for Plastics at 70°C.
RoHS 2 –Directive 2011/65/EU & Dangerous Substances – 76/769/EC
EU REACH Regulations – SVHC Current Candidate List
Bisphenol A – Absent
Latex – Absent
ADCF – Animal Derived Components Free (BSE/TSE). GMO Free
Endotoxins – Free
POP, PFOS, PFOA, APFO, DEHP, Phthalates, VOC Free
ODC – Class 1
Dodd-Frank Act – Gold, Coltan, Cassiterite, Wolframite and their derivatives originating from the Democratic Republic of Congo - Absent

It is the responsibility of the user to ensure that the products are suitable for their intended use

