


## Axygen® PCR 35 µm Aluminum Sealing Film For Thermal Sealing and General Sealing Purposes, Nonsterile

Product Number PCR-AS-200



Axygen® microplate sealing films are designed for applications ranging from PCR and Real-Time PCR to ELISA and cell culture. Axxygen's sealing films are the easiest and most cost-effective method for sealing plates. They can be used in wide temperature ranges, making them suitable for many applications.

- Provides exceptional seal on multiple well plates with a wide range of reagents and temperatures
- When used in conjunction with Axxygen's compression mats, seals during thermal cycling
- Ideal for light-sensitive samples
- Uniformly applied adhesive provides optimum well sealing

Product Number	PCR-AS-200
Qty./Pk	100 / Pk
Thickness	35 µm
Brand	Axygen®
Plate Format	Sealing films
Sterile	No