

TECHNICAL DATA SHEET

Part Number	11090-25G-HON
Description	Arsenazo III

Technical specifications :

- ✓ arsenazo iii, for spectrophotometric det. of th, zr, u, cd, zn, ca 25g
- ✓ CAS N°1668-00-4
- ✓ Molecular weight 776.37 g/mol
- ✓ Molecular formula C₂₂H₁₈As₂N₄O₁₄S₂
- ✓ no data available
- ✓ Flash point Not applicable
- ✓ Material form solid

- ✓ Storage condition Room Temperature
- ✓ Package type Glass Bottle
- ✓ Package size 25g

- ✓ GHS pictograms GHS09 ; GHS06
- ✓ H Statements H410 ; H301 ; H331
- ✓ P Statements P284 ; P280 ; P280 ; P260

Sales unit :

25g Arsenazo III
