

**1. Part Number:** 8111P-350

**2. Photograph:**



### 3. Product Information:

Description:	SteriWare® Powder Funnel - 350mm
Use:	Single use, disposable funnel
Material of Construction:	Polypropylene (PP) - virgin <ul style="list-style-type: none"> <li>• Conforms to FDA 21 CFR 177.1520</li> <li>• Conforms to EU Regulations 10/2011</li> <li>• Conforms to EC Regulations 1935/2004</li> <li>• BSE/TSE free</li> </ul>
Method of Construction:	Single piece, injection moulded
Maximum Diameter:	343mm
Overall Height:	264mm
Outside Diameter of Tip of Spout:	52mm
Inside Diameter of Tip of Spout:	48mm
Length of Spout:	43mm
Nominal Weight:	348g
Moulding Environment:	General Factory
Primary Packing Environment:	Class 10,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE bag)
Method of Sterilisation:	N/A
Number of Funnels per Box:	8

Recommended Storage Conditions:

Dry and ambient temperature

Shelf Life:

5 years from the date of manufacture

ATEX:

This product has not been ATEX tested

#### 4. Change Control:

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to

[change.control@sampling.com](mailto:change.control@sampling.com)

#### 5. Country of Origin:

Made and packed in the UK.

#### 6. Repeated Use:

This is a single-use, disposable product.

#### 7. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 23-Nov-22