

1. Part Number: M1013B-170

2. Photograph:



3. Product Information:

Description:

Metal Detectable Micro Spatula – 170mm - Blue

Material of Construction:

Polypropylene (PP)

- Metal detectable grade
- Only virgin material is used
- PP conforms to FDA 21 CFR 177.1520
- Blue colourant conforms to FDA 21 CFR 178.3297
- Filler conforms to FDA 21 CFR 186.1374
- Conforms to EU Regulations 10/2011
- Conforms to EC Regulations 1935/2004
- BSE/TSE free

Method of Construction:

Injection moulded as a single piece

Spatula Length:

175mm

Spatula Blade Length:

40mm

Blade Width:

13mm

Nominal Weight:

6g

Moulding Environment:

Class 100,000 Cleanroom

Packing Environment:

Class 10,000 Cleanroom

Individually Bagged?

Yes (heat sealed PE/PET bag)

Number of samplers per Pack:

4

Recommended Storage Conditions:

Dry and ambient temperature

Shelf Life:

5 years from the date of manufacture

ATEX:

This device has not been ATEX tested

4. Country of Origin:

Made and packed in the UK.

5. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 24-Jan-23

ISO
9001 : 2015
REGISTERED
Cert No. 14128455**Sampling Systems Ltd.**4-5 Forge Mills Park, Coleshill. Warwickshire. B46 1JH. UK.
Tel: +44(0)1675 466992 Email: info@sampling.com Web: www.sampling.com