



1. Part Number: 5065N-190

2. Photo:



Note: Telescopic Rod is sold separately.

3. Product Information:

Description	Deep Catch Net
Nominal Mesh Size:	2mm
Materials:	Nylon, Plastic
Max. Width:	500mm
Net Width:	440mm
Net Length:	190mm
Net Depth:	340mm
Nominal Weight:	330g
Recommended Storage Conditions:	Dry and ambient temperature

Note: Suitable for use with 5065A range of telescopic rods.

4. Approval:

Signed:

Name: Kuldeep Flora
Position: Quality Manager
Date: 5-Aug-22



ISO
9001 : 2015
REGISTERED
Cert No. 14128455

Sampling Systems Ltd.

4-5 Forge Mills Park, Coleshill, Warwickshire, B46 1JH, UK.
Tel: +44(0)1675 466992 Email: info@sampling.com Web: www.sampling.com