

**1. Part Number:** A324-120

**2. Photo:**



**3. Product Information:**

Description	Ladle 120ml
Nominal Volume:	120ml
Material of Construction:	204 stainless steel
Method of Construction:	Crevice free, all welds ground and polished
Surface Finish:	Better than 0.5 microns Ra
Handle Length:	270mm
Bowl Diameter:	82mm
Nominal Weight:	128g
Recommended Storage Conditions:	Dry and ambient temperature

Note: This product should be thoroughly cleaned before it is used for the first time

**4. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 3-Sep-21