

# Polyethylene glycol 6000

Ph. Eur., USP

(PEG 6000, Macrogol 6000)

CAS [25322-68-3]

Polyethylene glycol for chromatography, histology, microscopy and for special biochemical purposes.

Degree of polymerization ca. 140 - 170.

Average $M_r$	5400 - 6600
---------------	-------------

Hydroxyl value	17.0 - 21.0
----------------	-------------

EINECS: 500-038-2(NLP) • WGK: 1L • HS: 34042000

Storage Temperature: +15 °C to +30 °C