

# Polyethylene glycol 40 000

Ph. Eur., USP

(PEG 35000, Macrogol 35000)

CAS [25322-68-3]

Degree of polymerization ca. 800 - 900.

Polyethylene glycol for chromatography, histology, microscopy and for special biochemical purposes.

Average $M_r$	25 000 - 40 000
Heavy metals (Pb)	max. 5 ppm

EINECS: 500-038-2(NLP) • WGK: 1L • HS: 39072911

Storage Temperature: +15 °C to +30 °C