

Trichloroacetic acid, 20 % solution

(TCA, TCA solution)

Trichloroacetic acid is used for protein precipitation from biological samples and for fixation of native and IEF PAGE gels.

Aqueous solution. TCA (cat. no. 36910): 200 g/L.

DANGER

Hazard

Statements

Precautions

Reaction

H314 - H335 - H336 - H412

P273 - P280

P301 +P330 +P331 - P303 +P361 +P353 - P304 +P340 - P305
+P351 +P338

GGVSE/ADR: 8 III UN2564 • IATA: 8 III UN2564 • WGK: 2S • HS: 38229000

Storage Temperature: +15 °C to +30 °C