

# Technical Datasheet

---

**Article no.: 1197**

**2-Propanol puriss. (min. 99.8 %)**

CH<sub>3</sub>CH(OH)CH<sub>3</sub>

For laboratory use.

Parameter	Value
Appearance of solution	complies
Acidity	max. 0.001 %
Water (KF)	max.0.1 %
Refractive index (20 °C)	1.3757 – 1.3787
Residue on evaporation	max. 0.01 g/l
Colour (Apha)	max. 5

Version No. 02/2026-02-03