

Cryo-Apron®

Protection from Splashes and Contact Cold

Features

- Cryogenic protection for ultra-cold applications down to -196°C (-320°F)
- Made from state-of-the-art materials that are lightweight for comfort and ease of movement
- Multiple lengths available for maximum safety
- Adjustable straps at the neck and waist to ensure a proper fit
- Quick-release buckles at the neck and waist for swift on/off function
- Safety certified to meet EU standards: EN 511, EN 388, EN ISO 13688:2013
- **Applications:** dispensing/transferring cryogenic liquids, handling of cryogenic samples into or removing from storage vessels



APRON LENGTH				
	36" (91.4cm)	42" (106.7cm)	48" (121.9cm)	54" (137.2cm)
Apron	Blue ● CA36 Pink ● P-CA36	Blue ● CA42 Pink ● P-CA42	Blue ● CA48 Pink ● P-CA48	Blue ● CA54 Pink ● P-CA54



CE Cat. III
Pink Apron C E Pending