

Acu•Punch®

Sterile, Disposable Skin Biopsy Punch

The ONLY complete range of
14 sizes: from 1mm to 12mm.

The ONLY disposable skin biopsy
punch **made in the U.S.A.**

Made in U.S.A.

Now 14 Sizes



**1mm to
12 mm Sizes**



Acuderm® inc.
Real Service/Real Value



Intertek

CMDCAS




1-800-327-0015

Acu•Punch®

Now Available In 14 Sizes

1 mm 1.5 mm 2 mm 2.5 mm 3 mm 3.5 mm 4 mm 4.5 mm 5 mm 6 mm 7 mm 8 mm 10 mm 12mm



The ONLY Complete Range of Sizes

Physicians can choose a punch that is precisely sized to the area being biopsied or excised. Whether it's a skin biopsy, hair transplant or acne scar graft, Acu•Punch® allows the physician to remove exactly the amount of tissue desired.

Always Sharp, Always Sterile

Acu•Punch® is always razor sharp and ready to use. Each Acu•Punch® comes individually packaged in a sterile peel-apart pouch.

Ergonomic Design

Acu•Punch® handle is ribbed for maximum control, especially when wearing surgical gloves. The handle is hollow and open at both ends for easy removal of excised tissue for biopsy. Each punch size is identified on the handle.

Try a sampler pack

Choose our sampler pack (catalog number PSAM) and receive one each of all 14 sizes. Or choose our combination pack (catalog number PCOM) and receive 5 each of our most popular sizes – 2 mm, 3 mm, 3.5 mm, 4 mm, and 5 mm.

Only From Acuderm

Look for the white handle and the Acuderm name to make sure you're getting the genuine, made-in-U.S.A. Acu•Punch® from Acuderm. No other brand we know of is made to our exacting quality controls and standards.



Stock Numbers

Size	Box of 25	Box of 50
1mm	P125	P150
1.5mm	P1525	P1550
2mm	P225	P250
2.5mm	P2525	P2550
3mm	P325	P350
3.5mm	P3525	P3550
4mm	P425	P450
4.5mm	P4525	P4550
5mm	P525	P550
6mm	P625	P650
7mm	P725	P750
8mm	P825	P850
10mm	P1025	P1050
12mm	P1225	P1250
Mixed Sizes (25 total)		
5 each of 2, 3, 3.5, 4, 5mm	PCOM	Box of 25
Sampler (14 total)		
1 each size, incl. wrapped	PSAM14	Box of 14
Individual (min. order 5)		
Any size	PIND	N/A
Mixed Large Sizes (15 total)		
5 each of 4, 5, 6mm size	VETMIX	Box of 15



5370 NW 35 Terrace, Suite 106

Fort Lauderdale, FL 33309

Toll-Free: (800) 327-0015

Tel: (954) 733-6935 • Fax: (954) 486-3602

E-mail: cust-service@acuderm.com

Web: www.acuderm.com

Acuderm products meet applicable world-class standards, including

Quality Certified and Registered ISO 13485:2003

Certified to Medical Directive CE

CMDCAS, UL 60601-1

Distributed By: