

**Manufacturer** Paul Marienfeld GmbH & Co. KG  
 Am Wöllerspfad 4  
 D-97922 Lauda-Königshofen  
 Germany

Single Registration Number (SRN) DE-MF-000023724

**Product** See appendix list of products

Intended purpose Cover glasses for counting chambers and hemacytometers are accessories for in vitro diagnostic examinations in optical microscopy and are used to cover cell suspension in counting chambers and to limit the volume of cell suspension over the counting grid.

EU Declaration of conformity according to Regulation (EU) 2017/746 on in vitro diagnostic medical devices

Classification Class [A] according to annex VIII, [rule 5a] of regulation (EU) 2017/746

**Conformity Assessment Procedure** according to article 48(10):  
 Annex I, II and III of regulation (EU) 2017/746

We declare under our sole responsibility that the medical device covered by this declaration is in conformity with regulation (EU) 2017/746.

This declaration of conformity is valid until the issue of a new version, but not later than 2026-12-31 and includes all products manufactured in this period.

Lauda-Königshofen, 2025-06-30

Responsible Person acc. Article 15,  
 IVDR

CEO




Gerd Marienfeld

Harry Marienfeld



Marienfeld \* Am Wöllerspfad 4 \* 97922 Lauda-Königshofen \* Germany

<b>List of products</b>	
<b>Art. No.</b>	<b>Basic- UDI-DI</b>
03190000914	4250317304HGMSS3
03190000925	4250317304HGMSS3
03190000926	4250317304HGMSS3
03190000929	4250317304HGMSS3
0319000098	4250317304HGMSS3
0319000099	4250317304HGMSS3
0350000	4250317304HGMSS3
03500031	4250317304HGMSS3
0350009991	4250317304HGMSS3
0351000	4250317304HGMSS3
0352000	4250317304HGMSS3
03520011	4250317304HGMSS3
0357004	4250317304HGMSS3
0358004	4250317304HGMSS3