


**NEW**

# Thorbi® ID Plate



## Sample Identification for Cryoboxes

Insert plate is made of transparent polypropylene (PP) with alphanumeric or numeric ID code. Stable and just 0.5 mm thick. Simply place under the grid divider in cryoboxes, whether plastic or cardboard. Polypropylene is resistant to low temperatures down to -196 °C.

Pack contains 10 pieces in a bag with press closure.

### For Box 133 x 133 mm

Dimensions Thorbi® ID Plate: 125 x 125 x 0.5 mm

	ID Code	Cat. No.
<b>Alphanumeric</b>		
9x9	A1-I9	LCBC0909AN
10x10	A1-J10	LCBC1010AN
<b>Numeric</b>		
9x9	1-81	LCBC0909NU
10x10	1-100	LCBC1010NU












### For Box 136 x 136 mm

128 x 128 x 0.5 mm

	ID Code	Cat. No.
<b>Alphanumeric</b>		
	A1-I9	LCAC0909AN
	A1-J10	LCAC1010AN
<b>Numeric</b>		
	1-81	LCAC0909NU
	1-100	LCAC1010NU

Special designs are possible.





Further highlights and products can be found in our catalogues and at [www.Kryobox.de](http://www.Kryobox.de)



Please ask for a personal fitting option.

[www.Kryobox.de](http://www.Kryobox.de)