




clinisafe⁺
SHARPSAFE
5^{ÈME} GÉNÉRATION 

PRODUCT DATA SHEET

General data	
Item code	SHA-51912613
GTIN code	5056271508450
Product Name	Clinisafe® 5 th GEN. 8 liters
Description	Sharps container 8 liters composed by a body and a plain lid
Brand	Clinisafe®
Business Unit	Safety Solutions
Category of safety packaging	Sharps container for residual liquids
Indication	Segregation, protection, transport, disposal of solid sharps biohazard waste
Characteristics data	
Nominal volume	8 l
Total volume	8,64 l
Usefull volume	6,00 l
Maximum weight	2 Kg
Body background aspect	Neutral dark grey
Lid color	Yellow
Label or printing color	Yellow
Inner bag color	NC
Recognizable dominant color	Yellow
Marking languages	FR/NL
Dimensions (assembled) in mm :	
Long	315 mm
Wide	230 mm
Height	250 mm
Opening	270 x 170 mm
Fill line height	177 mm
Weight	535 g
Material	Polypropylene recycled 100% (body) + PP virgin 100% (lid)
Phtalates & Latex free	Yes
Carbon footprint reduction	NA
Autoclavable	Non autoclavable
Autoclavability conditions	-
Temporary closure	Yes
Permanent closure	Yes
Handles	1
Availability Accessories	Stands, trolley brackets (see brochure and compatibility tables)
	
Logistic data	
Carton quantity	20
Carton weight	10,75 Kg
Carton dimensions (LxWxH)	457 x 325 x 724 mm
Carton volume	107,5321 dm3
Palet quantity	240 units
Pallet weight	144 Kg
Pallet dimension (LxWxH)	120 x 80 x 230 cm
Certifications data	
ADR (UN3291)	Yes
EN ISO23907-1:2019	Yes
Quality mark	BSI Kitemark
TRBA250	Yes
NFX30-507:2018	NC
CE MARK	94/62/EC packaging and packaging waste
Quality management system :	
ISO 9001	Yes
ISO 13485	Yes
ISO 14001	Yes
ISO 45001	Yes
Manufacture data	
Made in	United Kingdom
By	Frontier Plastics Ltd
Marketed in	A/B/CH/D/F/LUX/NL
By	Hospidex Group