

Latex finger cots antistatic

- Powder-free
- Smooth surface
- ESD capability
- Recommendation: ISO 6



Article Data

Dutscher code	Item no.	Size	Diameter
165087	380015	S	15 mm
165088	380016	M	18 mm
165089	380017	L	21 mm
165090	380018	XL	24 mm

Technical Details

Specifications	
Material	100% Natural latex
Colour	Black
Thickness	Ca. 0.1 ± 0.02 mm
Length	Ca. 63 ± 2 mm
Surface resistivity	10 ⁶ Ω ~ 10 ⁸ Ω/square
Packaging	1,440 pieces per antistatic bag

Distribution

All statements without guarantee. They were compiled on the basis of extensive research. A legal liability can not be derived from this.