

1000 µL Entry/Exit Conductive Filter Tip Boxes

Specification Sheet (Part/REF # 235266)



Product Specifications

Platform Compatibility	Microlab VANTAGE Microlab STAR Microlab NIMBUS
Tip Type	Standard
Tip Volume	1000 µL
Maximum Tip Capacity	1000 µL
Filter	Filter
Sterile	Non-Sterile
Tip Color (Conductivity)	Black (Conductive)
Barcode Label Color	White
Tip Material	Polypropylene - PP
Pipetting Compatibility	1000 µL channel - 96 MultiProbe Head
Rack Material	Polypropylene - PP
Tip Recycling Code	5
Rack Recycling Code	5
Purity Level	Free of DNase, RNase, Pyrogen (endotoxins), Human DNA, and PCR inhibitors
Tips / Rack	96 Tips/Rack
SLAS ANSI	Yes
Tip Length	95.1 mm
Orifice Size	0.7 mm
Case Size	Case of 4,800 tips (5 tip boxes/stack; 10 E/E stacks per case)
Packaging Option	Entry/Exit CO-RE Tip Box Stacks

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/automated-liquid-handling/disposable-tips/1000-µl-entry-exit-conductive-filter-tip-boxes>