


PDS No. 1882xx/ 188271-TRI	PRODUCT DATA SHEET				Page 1 of 1
Revision 15	Tube, 15 ml, PP, Screw Cap, Sterile				
	Greiner Item-No. 1882xx /188271-TRI				
Valid for Item-No.:	188261 (sterile)	188271 (sterile)	188272 (sterile)	188280 (sterile)	188281 (sterile)
	188283 (sterile)	188285 (sterile)	188271-TRI (sterile)		

1.	Description / Specification	
1.1	Description	PP tube, 15 ml, with screw cap, conical bottom, graduation and writing area, sterile 188271-TRI: triple packed
1.2	Dimensions	See Customer Drawing
1.3	Volume	Working volume: 14 ml
1.4	Material / Resin	Tube: PP (Polypropylene), free of heavy metal Cap: HDPE (High Density Polyethylene), free of heavy metal
1.5	Colour	Tube: 188261 -271, -272, -271-TRI, -281, -285: natural ; print: graduation blue, writing area white 188280, -283: brown ; print: graduation blue, writing area white Cap: 188261, -271, -272, -271-TRI, -280, -283: blue 188281, -285: white
1.6	Sterilization	188261, -271, -280, -281, -283, -285: SAL10 ⁻³ 188272, -271-TRI: SAL10 ⁻⁶
1.7	Quality Control	- <u>Raw Material-Control</u> : physical testing - <u>Product-Control</u> : testing of attributive and variable characteristics in accordance with the valid specification
1.8	Intended Use	General laboratory product for the processing, centrifugation and storage of samples to be used by qualified personnel in a laboratory environment.
1.9	Other Information	For single use only

2.	Features	
2.1	Basic features	Free of detectable DNase/RNase, human DNA and pyrogens. Contents non-cytotoxic
2.2	Temperature range	For application: -80°C to +100°C
2.3	Autoclavability	Yes (without cap)
2.4	Centrifugation, max. RCF	4.000 x g: swinging-bucket rotor 15.000 x g: fixed-angle rotor
2.5	Chemical Resistance	See homepage: https://www.gbo.com/en_INT/know-how-services/download-center.html
2.6	Shelf life	5 years .
2.7	Other Information	-

3.	Packaging	188261, -281, -283	188271, -272, -280	188285	188271-TRI
3.1	Pieces / bag or rack	50 in rack / bag	100 / bag	100 / bag	5 / bag (triple packed)
3.2	Pieces / box	500 pieces (10 racks)	1.000	700	500
3.3	Lot-No.	E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)			
3.4	Other Information	Certificate of Quality to download			

4.	Other Information
4.1	-

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 14	Date 23 January 2025	Date 24 January 2025	Date 24 January 2025	
Date 27.10.2021	Name S. Kaelberer	Name Dr. T. Schreiber	Name Dr. C.-K. Chai	

DISCLAIMER: The description of a certain product can only be considered as a guidance because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.

INTERNAL