



TWIRL'EM®

TWIRL'EM®

The **TWIRL'EM®** provides a safe, contaminant-free and flexible container that ensures reliable test results. Its integrated secure closure system makes it a sealed and airtight solution. Proven efficient and versatile, it is used by many professionals in various industries.

The simple and economical tool that will support you in all your sampling activities.

Features and benefits

- Made of FDA compliant material
- Can accommodate solid, semi-solid and liquid samples
- Maximum bag opening to facilitate insertion of sample
- Various sizes and closure combinations available depending on the type of sample
- Suitable for use in the homogenizer
- Guaranteed sterility
- RNase, DNase and pyrogen free




Inspired by your challenges

LABPLAS

VISIT LABPLAS.COM FOR ADDITIONAL INFORMATION ABOUT OUR PRODUCTS | Tel: 1.450.649.7343



TWIRL'EM® BAG INSTRUCTION



Inspired by your challenges

LABPLAS

VISIT LABPLAS.COM FOR ADDITIONAL INFORMATION ABOUT OUR PRODUCTS | Tel: 1.450.649.7343

